



Microfiber Isolation Gown

2021 v.1

The Facts:

- Fluid repellent
- Fabric meets level 2
- Wrap-around design (available)
- Comfortable and lightweight
- Comfortable quick dry, stretch knit cuffs safely cover the wrist area



STYLE	FULL SWEEP	SLEEVE INSEAM	WEIGHT	DESIGN
EHX4200/WA	63"	24"	280 gr	Wrap-Around Coverage
EHX4200-WAT	63"	24"	280 gr	Wrap-Around Coverage with Traceability Grid
EHX4200-T	63"	24"	280 gr	Back-Tie with traceability Grid
GOWN/PPE-3XL	74"	24"	307 gr	Non reversible/Back-Tie.. Grey Binding



Contact your sales consultant

mipinc.com



MIP

The Softer Side of Healthcare