

INTEGRIDERM



SPECIFICATIONS

STYLE #
MIP6080

DIMENSIONS
35" x 80" x 5.5"

WARRANTY
5 Years

NOTE: Available in 78", 80", or 84" standard mattress lengths.



PRESSURE REDISTRIBUTION MATTRESS

| FEATURES | BENEFITS |
|--|--|
| Cross Die-Cut Topper | <ul style="list-style-type: none"> - Improves air flow - Reduces heat build up - Minimizes sheer during ingress and degress - Cross die-cut design - Exclusive to MIP |
| Full width topper | <ul style="list-style-type: none"> - Topper is full width of mattress for maximum sleeping surface - Reduces discomfort when sleeping on or near edge of mattress |
| Firm Base Layer of foam | <ul style="list-style-type: none"> - Provides increased support - Prolongs life of the mattress |
| Bolstered rails to the top of base layer | <ul style="list-style-type: none"> - Increases edge support - Decreases chances of entrapment - Provides more stability around the perimeter |
| Flex points in siderails | <ul style="list-style-type: none"> - Allows for smoother gatching - Siderail foam is not stressed and will not show wear (wrinkling) - Decreases risk of separation from center core |
| Fire Barrier Sleeve | <ul style="list-style-type: none"> - Adding the Fire Barrier sleeve to the foam mattress will meet more stringent fire standards, ie.TB603, CAL129, BFD IX-11 |
| 5 Year Full Warranty | <ul style="list-style-type: none"> - Great durability for maximum cost-savings - Cover and Foam warranted for 5 years <p>*Please note: This mattress does not include handles on the underside of cover</p> |

TEXTILE SOLUTIONS FOR YOUR EVERY DAY NEEDS