

INTEGRIDERM



SPECIFICATIONS

STYLE #

MIP9000/VBR

DIMENSIONS

35" x 80" x 7"

WARRANTY

7 Years

NOTE: Available in 78", 80", or 84" standard mattress lengths.



**PRESSURE REDISTRIBUTION
MATTRESS**

MIP9000/VBR

FEATURES	BENEFITS
Cross Die-Cut Gel Infused Visco Topper	<ul style="list-style-type: none"> - Enhances patient comfort - Improves air flow - Reduces heat build up - Minimizes sheer during ingress and degress - Cross die-cut design - Exclusive to MIP
Full width topper	<ul style="list-style-type: none"> - Topper is full width of mattress for added sleeping surface - Reduces discomfort when sleeping on or near edge of mattress
Cross Die-Cut Support Layer	<ul style="list-style-type: none"> - Further improves air flow
Higher grade of foam throughout the mattress	<ul style="list-style-type: none"> - Provides increased support - Prolongs life of the mattress
Bolstered rails to the top of the third layer	<ul style="list-style-type: none"> - Increases edge support - Decreases chances of entrapment - Provides more stability around the perimeter
Flex points in siderails	<ul style="list-style-type: none"> - Allows for smoother gatching - Siderail foam is not stressed and will not show wear (wrinkling) - Decreases risk of separation from center core
Fire Barrier Sleeve	<ul style="list-style-type: none"> - Adding the Fire Barrier sleeve to the foam mattress will meet more stringent fire standards, ie.TB603, CAL129, BFD IX-11
7 Year Full Warranty	<ul style="list-style-type: none"> - Excellent durability for maximum cost-savings - Cover and Foam warranted for 7 years

TEXTILE SOLUTIONS FOR YOUR EVERY DAY NEEDS