



Swift® UltraSlide® System

THE AFFORDABLE,
WASHABLE, IN-SITU,
BED REPOSITIONING
SYSTEM

2022v. 3

The Swift® UltraSlide® features 2 elements, a Swift® Slider & UltraSlide® bottom sheet, which work together to support repositioning and turning manoeuvres in bed.

A practical solution for care givers


- Saves time & energy (remains on the bed when not in use)
- Can be washed alongside regular bed linen & in temperatures up to 167°F (75°C)
- Fits most healthcare mattresses including standard & batriatric sizes, profiling & alternating air pressure reducing

Maximum patient care

- Soft touch, comfortable fabrics against the skin
- Breathable taffeta to reduce heat buildup
- Discreet design, looks like standard bedding when not in use



Swift®



REWASH	REUSE	REPEAT	REDUCE WASTE

Below the icons is a large green rectangular area.

At the bottom of the green area are four social media icons: Facebook (f), Twitter (bird), LinkedIn (in), and YouTube (You Tube).

MIP

The Softer Side of Healthcare