



# SYNTHETIC UNDERPAD

2021 v.1

## The Facts

- Bonded quilted synthetic underpad with lightly brushed barrier
- 100% polyester, napped white Sonoma surface
- Waved quilted surface provides increased comfort and improved durability
- Needled non-woven soaker
- High performance at value price



## Specifications

STYLE #	DIMENSION	WEIGHT	BARRIER	SOAKER
11026/LV	34" x 36"	15.3 oz	Dark Blue	7 oz/y <sup>2</sup>



Contact your sales consultant

[mipinc.com](http://mipinc.com)



# MIP

The Softer Side of Healthcare